

Overview of NMR-STAR Objects

Entry

Data:

- BMRB ID : 15000
- Saveframes: [<bmr.b.saveframe 'entry_information'>, <bmr.b.saveframe 'citation_1'>, ...]

Methods:

- Construct:
 - From scratch, from a file, from a string, from the public database of entries
- Add and remove saveframes
- Search for saveframes by name or by category
- Generate a string of the entry in NMR-STAR format
- Validate and compare to other entries
- Print structure of entry as tree

Saveframe

Data:

- Name : 'entry_information'
- Tag prefix: '_Entry'
- Tags: {'Sf_category':'entry_information', 'Sf_framecode':'entry_information', ...}
- Loops: [<bmr.b.loop '_Entry_author'>, <bmr.b.loop '_SG_project'>, ...]

Methods:

- Construct:
 - From scratch, from a file, from a string
- Add and remove tags
- Add and remove loops
- Search for loops by category
- Generate a string of the saveframe in NMR-STAR format
- Generate a CSV containing the data in the tags
- Validate and compare to other saveframes
- Print tree structure of saveframe

Loop

Data:

- Category : '_Entry_author'
- Columns:
 - ['Ordinal', 'Given_name', 'Family_name', 'First_initial', 'Middle_initials', 'Entry_ID1']
- Data:
 - [['1', 'Claudia', 'Cornilescu', '', 'C.', '15000'], ['2', 'Gabriel', 'Cornilescu', '', '', '15000'], ...]

Methods:

- Construct:
 - From scratch, from a file, from a string
- Add columns
- Add, remove, and renumber data
- Fetch data by column and/or fetch row matching certain data
- Generate a string of the loop in NMR-STAR format
- Generate a CSV containing the data
- Validate and compare to other loops

An example entry

<bmr.entry '15000'>

BMRB ID and saveframe list:

**<bmr.saveframe
'entry_information'>**

Tag dictionary, saveframe
name, tag prefix, and loop list:

<bmr.loop '_Entry_author'>

<bmr.loop '_SG_project'>

...

**<bmr.saveframe
'citation_1'>**

Tag dictionary, saveframe
name, tag prefix, and loop list:

<bmr.loop '_Citation_author'>

**<bmr.saveframe
'assembly'>**

Tag dictionary, saveframe
name, tag prefix, and loop list:

<bmr.loop '_Entity_assembly'>

**<bmr.saveframe
'F5-Phe-cVHP'>**

Tag dictionary, saveframe
name, tag prefix, and loop list:

<bmr.loop '_Entity_db_link'>

<bmr.loop '_Entity_comp_index'>

<bmr.loop '_Entity_poly_seq'>

...

More detailed example of a saveframe:

<bmrbl.saveframe 'assigned_chem_shift_list_1'>

Tag_prefix: '_Assigned_chem_shift_list'

Tags: {'Sf_category':'assigned_chemical_shifts',
'Sf_framecode':'assigned_chem_shift_list_1','Entry_ID':'15000','ID': '1'), ...}

Loops: [<bmrbl.loop '_Chem_shift_experiment'>, <bmrbl.loop
'_Atom_chem_shift'>]

<bmrbl.loop '_Chem_shift_experiment'>

Category: '_Chem_shift'

Columns: ['Experiment_ID', 'Experiment_name', ...]

Data: 2x2 array of data

In NMR-STAR format:

```
loop_  
  _Chem_shift_experiment.Experiment_ID  
  _Chem_shift_experiment.Experiment_name  
  _Chem_shift_experiment.Sample_ID  
  _Chem_shift_experiment.Sample_label  
  _Chem_shift_experiment.Sample_state  
  _Chem_shift_experiment.Entry_ID  
  _Chem_shift_experiment.Assigned_chem_shift_list_ID  
  
  . '2D 1H-15N HSQC' 1 $unlabeled_sample isotropic 15000 1  
  . '2D 1H-13C HSQC' 1 $unlabeled_sample isotropic 15000 1  
  . '2D COSY' 1 $unlabeled_sample isotropic 15000 1  
  . '2D TOCSY' 1 $unlabeled_sample isotropic 15000 1  
stop_
```

<bmrbl.loop '_Atom_chem_shift'>

Some sample code:

Start the python interpreter and load the module

```
$: python
>>> import bmrB
>>> ent15000 = bmrB.entry.fromDatabase(15000)
>>> ent15000.printTree()
<bmrB.entry '15000'>
  [0] <bmrB.saveframe 'citation_1'>
    [0] <bmrB.loop '_Citation_author'>
  [1] <bmrB.saveframe 'assembly'>
    [0] <bmrB.loop '_Entity_assembly'>
  [2] <bmrB.saveframe 'F5-Phe-cVHP'>
    [0] <bmrB.loop '_Entity_db_link'>
    [1] <bmrB.loop '_Entity_comp_index'>
    [2] <bmrB.loop '_Entity_poly_seq'>
  [3] <bmrB.saveframe 'natural_source'>
    [0] <bmrB.loop '_Entity_natural_src'>
  [4] <bmrB.saveframe 'experimental_source'>
    [0] <bmrB.loop '_Entity_experimental_src'>
  [5] ....
```

There is a shorthand way to access saveframes by name and loops by category.

```
>>> ent15000['entry_information']
<bmrB.saveframe 'entry_information'>
>>> ent15000['entry_information']['_Entry_author']
<bmrB.loop '_Entry_author'>
```

```
>>> ent15000_2 = bmrB.entry.fromDatabase(15000)
>>> del ent15000_2['entry_information']
>>> bmrB.diff(ent15000,ent15000_2)
```

The number of saveframes in the entries are not equal: 25 vs 24
No saveframe with name 'entry_information' in other entry.

Let's look at a loop's column headers and its data

```
>>> ent15000['entry_information']['_Entry_author'].columns
['Ordinal', 'Given_name', 'Family_name', 'First_initial',
'Middle_initials', 'Family_title', 'Entry_ID']
>>> ent15000['entry_information']
['_Entry_author'].getDataByTag('Given_name')
[['Claudia'], ['Gabriel'], ['Erik'], ['Samuel'], ['John']]
```

Write the entry to disk in NMR-STAR format

```
>>> with open('/tmp/entry','w') as tmp_file:
...     tmp_file.write(str(ent15000))
```

Get a list of validation errors

```
>>> ent15000.validate()
["Tag '_Entry.Type' not found in schema. Line 4.", "Value cannot be NULL
but is: '_Chem_comp.Provenance': '.' on line 4.", ...]
```

Here is how to create a loop from scratch

```
>>> new_loop = bmrbl.loop.fromScratch()
>>> new_loop.addColumn("_Example.day")
>>> new_loop.addColumn("month")
>>> new_loop.addColumn("year")
>>> new_loop.addColumn("event")
```

You could replace above with:

```
>>new_loop.addColumn(["_Example.day","month","year","event"])
```

Add data to the loop

```
>>> new_loop.addData([25,11,13,"Nothing exciting happened."])
```

Notice that data is automatically encapsulated as necessary to meet STAR format (quotes around the data containing a space). You never have to worry about encapsulating data you insert.

```
>>> print new_loop
```

```
loop_
  _Example.day
  _Example.month
  _Example.year
  _Example.event

    25    11    13    'Nothing exciting happened.'
stop_
```

Add the loop to the entry_information saveframe

```
>>> ent15000['entry_information'].addLoop(new_loop)
```

```
>>> ent15000['entry_information'].printTree()
```

```
<bmrbl.saveframe 'entry_information'>
  [0] <bmrbl.loop '_Entry_author'>
  [1] <bmrbl.loop '_SG_project'>
  [2] <bmrbl.loop '_Struct_keywords'>
  [3] <bmrbl.loop '_Data_set'>
  [4] <bmrbl.loop '_Datum'>
  [5] <bmrbl.loop '_Release'>
  [6] <bmrbl.loop '_Related_entries'>
  [7] <bmrbl.loop '_Example'>
```

```
>>> print ent15000['assembly'].tags
```

```
OrderedDict([('Sf_category', 'assembly'), ('Sf_framecode', 'assembly'),
('Entry_ID', '15000'), ('ID', '1'), ('Name', 'F5-Phe-cVHP'), ('BMRB_code',
 '.'), ('Number_of_components', '1'), ...])
```

```
>>> ent15000['assembly'].getTag('ID')
```

```
1
```

```
>>> ent15000['assembly'].printTree()
```

```
<bmrbl.saveframe 'assembly'>
  [0] <bmrbl.loop '_Entity_assembly'>
```

Get a saveframe in NMR-STAR format

```
>>> print ent15000['assembly']
```

```
save_assembly
```

| | |
|-------------------------------------------|-------------|
| _Assembly.Sf_category | assembly |
| _Assembly.Sf_framecode | assembly |
| _Assembly.Entry_ID | 15000 |
| _Assembly.ID | 1 |
| _Assembly.Name | F5-Phe-cVHP |
| _Assembly.BMRB_code | . |
| _Assembly.Number_of_components | 1 |
| _Assembly.Organic_ligands | . |
| _Assembly.Metal_ions | . |
| _Assembly.Non_standard_bonds | . |
| _Assembly.Ambiguous_conformational_states | . |
| _Assembly.Ambiguous_chem_comp_sites | . |
| _Assembly.Molecules_in_chemical_exchange | . |
| _Assembly.Paramagnetic | no |
| _Assembly.Thiol_state | 'all free' |
| _Assembly.Molecular_mass | . |
| _Assembly.Enzyme_commission_number | . |
| _Assembly.Details | . |
| _Assembly.DB_query_date | . |
| _Assembly.DB_query_revised_last_date | . |

```
loop_
```

| |
|--------------------------------------------------|
| _Entity_assembly.ID |
| _Entity_assembly.Entity_assembly_name |
| _Entity_assembly.Entity_ID |
| _Entity_assembly.Entity_label |
| _Entity_assembly.Asym_ID |
| _Entity_assembly.PDB_chain_ID |
| _Entity_assembly.Experimental_data_reported |
| _Entity_assembly.Physical_state |
| _Entity_assembly.Conformational_isomer |
| _Entity_assembly.Chemical_exchange_state |
| _Entity_assembly.Magnetic_equivalence_group_code |
| _Entity_assembly.Role |
| _Entity_assembly.Details |
| _Entity_assembly.Entry_ID |
| _Entity_assembly.Assembly_ID |

| | | | | | | | | | |
|----|---|-------------|-------|---------------|---|---|-----|--------|----|
| no | 1 | F5-Phe-cVHP | 1 | \$F5-Phe-cVHP | K | . | yes | native | no |
| | . | . | 15000 | 1 | | | | | |

```
stop_  
save_
```

Reading Spectral Peaks from a NMR-STAR file:

```
# First load the file and get a list of the peak list saveframes
```

```
>>> ent6577 = bmr.entry.fromDatabase(6577)
```

```
>>> spectral_peaks = ent6577.getSaveframesByCategory('spectral_peak_list')
```

```
# Lets look at how many spectral peak list saveframes we have
```

```
>>> spectral_peaks
```

```
[<bmr.saveframe 'peak_list_1'>, <bmr.saveframe 'peak_list_2'>,  
<bmr.saveframe 'peak_list_3'>]
```

```
# For this demo we'll just look at one individual peak list
```

```
>>> peak1 = spectral_peaks[0]
```

```
# We can see what loops this peak list saveframe contains
```

```
>>> peak1.printTree()
```

```
<bmr.saveframe 'peak_list_1'>
```

```
  [0] <bmr.loop '_Spectral_peak_software'>
```

```
  [1] <bmr.loop '_Peak_general_char'>
```

```
  [2] <bmr.loop '_Peak_char'>
```

```
  [3] <bmr.loop '_Assigned_peak_chem_shift'>
```

```
# Let's see what the _Peak_char loop looks like in NMR-STAR format
```

```
>>> print peak1['_Peak_char']
```

```
loop_
```

```
  _Peak_char.Peak_ID
```

```
  _Peak_char.Spectral_dim_ID
```

```
  _Peak_char.Chem_shift_val
```

```
  _Peak_char.Chem_shift_val_err
```

```
  _Peak_char.Line_width_val
```

```
  _Peak_char.Line_width_val_err
```

```
  _Peak_char.Phase_val
```

```
  _Peak_char.Phase_val_err
```

```
  _Peak_char.Decay_rate_val
```

```
  _Peak_char.Decay_rate_val_err
```

```
  _Peak_char.Coupling_pattern
```

```
  _Peak_char.Derivation_method_ID
```

```
  _Peak_char.Entry_ID
```

```
  _Peak_char.Spectral_peak_list_ID
```

| | | | | | | | | | | | | |
|---|---|-------|---|---|---|---|---|---|---|---|------|---|
| 1 | 1 | 9.857 | . | . | . | . | . | . | . | . | 6577 | 1 |
| 1 | 2 | 4.922 | . | . | . | . | . | . | . | . | 6577 | 1 |
| 2 | 1 | 9.857 | . | . | . | . | . | . | . | . | 6577 | 1 |
| 2 | 2 | 2.167 | . | . | . | . | . | . | . | . | 6577 | 1 |

```
...
```

```
stop_
```

Reading Spectral Peaks from a NMR-STAR file (cont.):

That is more information than we want right now. Lets get just the columns we need (we'll get a list of lists, each inner list corresponds to a row).

```
>>> our_data = peak1['_Peak_char'].getDataByTag(
    ['Peak_ID', 'Chem_shift_val'])
```

```
>>> print our_data
[['1', '9.857'],
 ['1', '4.922'],
 ['2', '9.857'],
 ['2', '2.167'],
 ['3', '9.855'],
 ...]
```

Excellent! Now we can iterate through each spectral peak and corresponding shift easily. The data is stored as a python list of lists (2 dimensional array) and we can modify or access it any of the normal ways python allows.

```
>>> for x in our_data:
...     print "Sprectral chemical shift value is: " + str(x[1])
...
Sprectral chemical shift value is: 9.857
Sprectral chemical shift value is: 4.922
Sprectral chemical shift value is: 9.857
...
```

It is also easy to dump the table in a loop as a CSV

```
>>> print peak1['_Peak_char'].getDataAsCSV()
Peak_ID,Spectral_dim_ID,Chem_shift_val,Chem_shift_val_err,Line_width_val,Line_width_val_err,Phase_val,Phase_val_err,Decay_rate_val,Decay_rate_val_err,Coupling_pattern,Bounding_box_upper_val,Bounding_box_lower_val,Bounding_box_range_val,Details,Derivation_method_ID,Entry_ID,Spectral_peak_list_ID
1,1,9.857,,,,,,,,,,,,,6577,1
1,2,4.922,,,,,,,,,,,,,6577,1
2,1,9.857,,,,,,,,,,,,,6577,1
2,2,2.167,,,,,,,,,,,,,6577,1
3,1,9.855,,,,,,,,,,,,,6577,1
3,2,1.994,,,,,,,,,,,,,6577,1
4,1,9.854,,,,,,,,,,,,,6577,1
```